

Intrinsic Optical Signal Imaging Code Instruction Manual

BMEG CAPSTONE 2024-2025

Client: Dr. Bridi, Rockefeller Neuroscience Institute

Dr. Anika Pirkey: [Anika Pirkey - arcoolbaugh@mix.wvu.edu](mailto:arcoolbaugh@mix.wvu.edu)

Jada Mullins

Hunter Jackson

Lauren Musbach

Riley Nett

Neopixel Code: NeopixelApp.mlapp

Instruction is made for Adafruit Neopixel LEDs with Integrated Drivers.
Adapted from the Zeiger Lab github linked in instructions.

Wire the Arduino Uno Rev 3

1. Connect an electrical wire (White) for input to the Neopixel to Digital Pin 6 of the Arduino Uno
 - a. Solder the wire end connected to the Neopixel ring to the base of the ring in points labelled "IN".
2. Connect a 2nd electrical wire (Black) for ground to the GND pin of the Arduino Uno.
 - a. Solder the end of the electrical wire to the point in the Neopixel ring labelled 'GND'.
3. Connect a 3rd electrical wire (Red) for the 5V connection input to the Digital pin '5V' of the Arduino Uno.
 - a. Solder the end of the electrical wire to the point in the Neopixel ring labelled 'PWR'.

Wire the Neopixel

4. Take pieces of pre insulated electrical wire (white, black, and red) that are of an approximate length of 30 cm.
5. Solder the other ends of each wire to the Neopixel Ring as follows:
 - a. White to "IN" hole
 - b. Black to the "GND" hole (there are two G holes, either is fine)
 - c. Red to the "PWR" hole (there are two V+ holes, either is fine)
 - d. Note: It is easiest to insert the wire from the front of the Neopixel (where the LEDs are located) and solder onto the back

Assemble the Neopixel Jewel Holder

6. Bolt down metal piece to the device where you wish for the gooseneck to begin.
7. Attach a gooseneck using the magnetic base to the metal piece. (Note: Aluminum will not work well for this, the metal needs to have a strong magnetic connection)
8. Cut a piece of double sided adhesive to fit the end portion of the gooseneck.
9. Use double sided adhesive to attach the Neopixel Jewel to the tip of the gooseneck.

Install Software

10. Connect a USB Type A/B cable between a computer and the Arduino Uno Rev3
11. Download and install the Arduino IDE (<https://www.arduino.cc/en/software>)
12. Download and install the Adafruit Neopixel Library (https://github.com/adafruit/Adafruit_NeoPixel)
13. Download and install MATLAB (<https://matlab.mathworks.com>)
14. Download and install the MATLAB Arduino Support Package (<https://www.mathworks.com/hardware-support/arduino-matlab.html>)
15. Download and install the MATLAB Neopixel Add-On Library for Arduino (<https://www.mathworks.com/matlabcentral/fileexchange/72707-neopixel-add-on-library-for-arduino>)


16. Clone our intrinsic signal imaging repository and add to your MATLAB path (<https://github.com/jadamullns/Bridi-Labs-IOSI>)
17. Determine which port your Arduino is connected to. You may need to change a line of code to ensure this is correct
 - a. Open “Device Manager” in Windows and expand the “Ports (COM & LPT)” list (<https://www.mathworks.com/help/supportpkg/arduinoio/ug/find-arduino-port-on-windows-mac-and-linux.html>)
 - i. Find the Arduino Uno and make note of the port listed (e.g. “COM5”)
 - b. Adjust the first input variable in line 31. For example:
 - i. `a = arduino('COM3', 'Uno', 'Libraries', 'Adafruit/NeoPixel');`
 - c. If you have adapted this protocol to use a different Neopixel Ring, you may need to further adjust the code to match the specifications of your Neopixel Ring
18. From the MATLAB window, open the app to control the Neopixel illumination.

Running Matlab App

19. Select Run icon on Matlab user interface to begin the app connection to Neopixel Ring
20. Select color of LED by selecting button on the app with the desired color labeled
21. Select Color intensity by typing in intensity value (from 1-100% Intensity)
22. Adjust throughout stimulation as indicated by the trial timing instruction guide.

IOSGUI Code: iosgui.mlapp

Adapted from the Zeiger Lab github linked in instructions.



Quickstart Guide

1. Open MATLAB and type “iosgui” to open the image processing application
2. Edit the fields in “Image Acquisition Settings” to match the settings used for acquiring IOSI images
3. Click “Import Imaging Files” and select a group of images (or an image stack) for processing
4. Click “Process Files” to start image processing and generate a scaled $\Delta R/R$ image
5. Click “Import Vasculature Image” to import an image of cortical vasculature taken contemporaneously with the IOSI images
6. *Optional* Click “Add Mask” and re-size the displayed ellipse to mask pixels outside the area of the ellipse
7. Click “Generate Z-score Overlay” or “Generate % Threshold Overlay” to generate an image with IOSI signals binarized and overlaid onto the cortical vasculature image


Installation

- 1) Clone the Intrinsic-Signal-Imaging repository from Github:
<https://github.com/zeigerlab/Intrinsic-Signal-Imaging>

- a. For detailed instructions, see <https://docs.github.com/en/repositories/creating-and-managing-repositories/cloning-a-repository?tool=desktop>
- 2) Add the folder containing the repository to your MATLAB path
 - a. For detailed instructions, see: https://www.mathworks.com/help/matlab/matlab_env/add-remove-or-reorder-folders-on-the-search-path.html
- 3) Install the “Image Processing Toolbox” in MATLAB:
<https://www.mathworks.com/products/image.html>

Image Acquisition Settings

- Number of Trials
 - o Total number of trials to analyze. A single trial consists of a baseline period, a stimulus (whisker deflection, drifting visual grating, etc.), and a post-stimulus period. Our IOSI experiment consists of 5 trials. If you have a different trial structure (for example, you only want to analyze every 3rd trial), prior to analysis you will need to create a folder containing that particular subset of images, starting at the first trial to be analyzed, ordered and numbered consecutively.
- Trial Length
 - o The total number of frames (images) acquired for each individual trial. Include all acquired frames, even if some may not contribute to the analysis. For example, if



you recorded for 4 seconds at 10 frames per second (or 10 Hz) for each trial, the trial length in frames will be 40. The ThorCam software takes images at 1 frame per second.

- Baseline Period
 - o The number of frames acquired before the stimulus of interest. For example, if the stimulus begins after 1 second of image acquisition at 10 Hz, the baseline period in frames will be 10. These images will be averaged on each trial to create a mean pre-stimulus reflectance image for the calculation of change in reflectance ($\Delta R/R$) values.
- Intrinsic Signal Period
 - o The number of frames (images) acquired during the period over which you would like to calculate intrinsic signals. This period is typically 1-1.5 seconds, but may

vary depending on the particular experiment and the specific intrinsic signal you are interested in quantifying. These images will be temporally binned (see below), normalized to the mean pre-stimulus image, and summed to calculate $\Delta R/R$.




Figure 2. Structure of an individual trial. Each trial is expected to consist of consecutively acquired frames comprising a pre-stimulus baseline period running from acquisition start to stimulus onset, an offset (optional) that allows for measurement of a desired phase of intrinsic signals, and the intrinsic signal period over which change in reflectance values will be calculated.

- Stimulus Offset
 - The number of frames (images) between the onset of the stimulus of interest and the beginning of the period over which you would like to calculate intrinsic signals.
- Temporal Binning
 - The number of frames (images) over which to bin consecutively acquired images from the intrinsic signal period. Temporal binning can be helpful to reduce noise, increasing signal-to-noise ratio. Generally, higher acquisition rates can benefit from temporal binning.
- Spatial Binning
 - A multiplicative factor by which to spatially downsample images using a bilinear interpolation. For example, a 1000 x 1000 pixel image with a spatial binning factor of 0.25 will become a 250 x 250 pixel image. Spatial downsampling can improve signal to noise ratio, but at the expense of spatial resolution. We recommend a spatial binning factor of 0.25.

File Selection

- Import Imaging Files
 - Opens a window for file selection. Currently, the application can accept grayscale tiff image files in one of two formats:
 1. Individual sequentially numbered files that share a common base name followed by an image number (such as “image_001.tif, image_002.tif, image_003.tif”).
 - a. Image numbers should span trials – if each trial is 40 frames long, the first frame of the second trial will be image 41
 - b. Hold “Control” or “Command” to select all images across all trials to be analyzed.
 2. A multi-page tiff in which all images across all trials are sequentially acquired and stored together in one file. This is the one we recommend.

- Process Files
 - This button will calculate intrinsic signals using the images selected previously and the settings defined in the input boxes in the “Image Acquisition Settings” section of the application.
 - A calculated $\Delta R/R$ image (scaled from minimum intensity to maximum intensity) will be displayed in the “Displayed Images” section of the application.
 - The scaled image and a MATLAB “.mat” file with raw data will be automatically saved in the directory containing the images selected for analysis. If scaled image and analysis MATLAB files already exist in the directory, a timestamp will be appended to the file names to prevent overwriting of data.
- Import Vasculature Image (Optional)
 - Opens a window for file selection. Choose a single grayscale tiff image of the cortical surface vasculature, typically acquired either immediately before or after an IOSI experiment.
 - The chosen vasculature image will be displayed in the “Displayed Images” section of the application.
 - In subsequent steps (see below) $\Delta R/R$ values can be binarized according to a threshold and overlaid onto this image to generate a map of intrinsic signals that can be localized according to the cortical vasculature.
- Add Mask (Optional)
 - This allows you to define areas of the acquired images to be excluded when generating overlay activity maps.
 - An ellipse will be displayed on the vasculature image in the “Displayed Images” section of the application.
 - Adjust and move the ellipse to fit the area of interest. The mask will automatically be updated each time the ellipse is adjusted or moved. Signals within the ellipse will be included and all signals outside of the ellipse will be excluded when generating overlay activity maps. This is useful, for example, for limiting signals to the area of a circular cranial window and excluding artifacts caused by surrounding hardware.

Program Messages

- Progress updates, error messages and warnings will be displayed in this area

Displayed Images

- Image Axes
 - Scaled $\Delta R/R$ images, imported vasculature images, and overlay images will be displayed on the axes in this area
- Adjust Contrast
 - This tool can be used to adjust the contrast or brightness of displayed images
 - When used with grayscale images (scaled $\Delta R/R$ images or imported vasculature images) a new window will open in which you can interactively adjust contrast limits. When ready, click the “Adjust Data” button.



Figure 3. The “Adjust Contrast” window.

- The adjust contrast tool is not compatible with RGB images (such as the overlay images generated by the app, see below). In this case, the image brightness will be increased iteratively by 5% each time the “Adjust Contrast” button is pressed.
- Adjustments are for display purposes only and are not automatically saved. If you would like to save an image with adjusted contrast, hover over the upper right corner of the Image Axes and use the image toolbar to save manually.
- Generate Z-score Overlay




Figure 4. The image axes toolbar. To manually save an image, hover over the upper right corner of the axes to access the toolbar.

- o This button will create an image of IOSI signals overlaid onto an image of the cortical vasculature using a Z-score threshold
 - Calculated $\Delta R/R$ values will be Z-scored and binarized according to the threshold set in the input box next to the button
 - In our experience Z-scores of -2 to -3 work well for first phase “initial dip” intrinsic signals. Smaller (more negative) Z-scores will exclude more signal, keeping just the strongest signals. Larger (more positive) Z-scores will include more signal and therefore potentially background noise. If you are measuring a later phase

intrinsic signal (such as the positive-going 2nd phase) these ranges will need to be adjusted accordingly.

- Signals outside of any mask previously set will be excluded
- The binarized signal will be pseudo-colored and overlaid onto a previously selected vasculature image. Both the binary image and the overlay will be automatically saved in the directory containing the images selected for analysis. If binary and overlay image files already exist in the directory, a timestamp will be appended to the file names to prevent overwriting of data.
- Generate % Threshold Overlay
 - o This button will create an image of IOSI signals overlaid onto an image of the cortical vasculature, using a threshold set according a percentage of the total signal range
 - A threshold, based on the full range of $\Delta R/R$ values (maximum minus minimum) will be calculated according to the percentile set in the input box next to the button. The $\Delta R/R$ values will then be binarized according to the threshold, keeping only the specified percentage of the overall range of values
 - In our experience percentiles of 30-60% work well for first phase “initial dip” intrinsic signals. Smaller percentiles will exclude more signals, keeping just the strongest signals. Larger percentiles will include more signal and therefore potentially background noise. If you are measuring a later phase intrinsic signal (such as the positive-going 2nd phase) these ranges will need to be adjusted accordingly.
 - Signals outside of any mask previously set will be excluded
 - The binarized signal will be pseudo-colored and overlaid onto a previously selected vasculature image. Both the binary image and the overlay will be automatically saved in the directory containing the images selected for analysis. If binary and overlay image files already exist in the directory, a time stamp will be appended to the file names to prevent overwriting of data.

Visual Stimulus shown on Tab 2